

NITROUSEAL® NITROUS ANALGESIA SYSTEM

The premium and only complete nitrous oxide system available.



Mouthpiece and Nitrouseal® circuit ideal for facial procedures



Nitrouseal® mask and circuit ideal for body procedures



Nasal hood ideal for peri-oral and neck procedures

Nitrouseal® Benefits

Benefits to Practitioners	Benefits to Patients
More comfortable patients: increase referrals, more likely to return for additional treatments, improve reputation of practice, easier to perform procedures, more revenue for practice	Pain and anxiety relief: non-opioid
Safe for patient: does not depress respiration or airway reflexes. Pt's awake and can communicate during procedure	Self-administered: easy to use
Safe for staff: all waste nitrous is scavenged and eliminated from procedure room so avoid risks of chronic exposure	Very rapid onset and offset
Easy to learn and easy to use: Nitrouseal® provides on-site training and is available to answer questions and provide additional on-line help	Patients can drive themselves home: nitrous eliminated so rapidly that patients are ready to go within 30 minutes
Minimal sedation: avoids regulations governing moderate sedation or deeper levels of anesthesia.	
Titrateable: to the individual patient to avoid the side effects of too much nitrous	

30-day money back guarantee if not completely satisfied (minus \$250 for shipping and any consumables)



5901 Sun Blvd., Ste. 200
St. Petersburg, FL 33715

www.sedationsystems.com | (888) 959-5288

One or more of the following U.S. Patents may be applicable to this product: 8,336,549; 8,826,905; D697,202; D658,284; D615,643; D598,537. International and U.S. patents pending.

©2020, Sedation Systems, LLC. All rights reserved.




Adjustable flowmeter from 0-70% nitrous (*can be limited to 50% by simple programming*)
100% oxygen at the end of procedure



Only FDA cleared portable scavenger available in the United States
(CDC, NIOSH, and OSHA all recommend scavenging to eliminate nitrous)

Only Complete Nitrous System Available



Patented mask and circuit tested and proven to minimize room contamination



Nitrouseal® Comparison	ProNox® or Nitronox™
BENEFITS	
1. Adjustable nitrous concentration	1. Locked at 50:50 only
2. Titrate to the individual; avoids unpleasant side effects of too much nitrous	2. No
3. Can give 100% of O2 at beginning and end of each procedure [very important safety feature]	3. No
4. Scavenging FDA cleared: eliminates room contamination with nitrous [very important for staff safety] – all exhaled air from patient scavenged [very important in COVID 19 era]	4. No
5. Patented mask and circuit designed by an anesthesiologist for <u>safety, comfort</u> , and to <u>minimize room contamination</u> [tested at two large pediatric hospitals and proven to minimize nitrous contamination]	5. Generic mask
SAFETY	
1. Alarms inform which gas has run out	1. Does not tell you which gas has run out
2. Four tank stand – always have a backup tank	2. Only has one of each with no back up
REGULATIONS	
In some states only Nitrouseal® can be used. For example, Delaware only allows nitrous below 50%. Nitrouseal® can be used because it can administer any concentration of nitrous between 0-49%.	